

Pe Rabbit Anti-Mouse Integrin Alpha 1/Cd49A Mab For Fc - Q3V3R4

Item Number: CM0008398

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse Integrin alpha 1/CD49a (VLA-1) for flow cytometry. Detects cell surface receptor for laminin and collagen adhesion signaling.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse Integrin alpha 1/CD49a mAb
Remarks/Alias	Vla1, CD49A, E130012M19Rik
Species	Mouse
GeneID (Human)	3672
GeneID	109700
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-1141 of mouse Integrin alpha 1/CD49a (NP_001028400.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q3V3R4
Protein Weight	130kDa
Shipping	Ice bag