


Pe Rabbit Anti-Mouse Lamp1/Cd107A Monoclonal Antibody For Fc - P11438

Item Number: CM0012103



Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse LAMP1/CD107a for flow cytometry. Validated for mouse samples. LAMP1 is a lysosomal membrane glycoprotein involved in lysosome biogenesis, autophagy, and NK cell cytotoxicity.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse LAMP1/CD107a mAb
Remarks/Alias	P2B; Perk; LGP-A; CD107a; Lamp-1; LGP-120
Species	Mouse
Gene ID (Human)	3916
Gene ID	16783
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P11438
Protein Weight	43kDa
Shipping	Ice bag