

Pe Rabbit Anti-Mouse Ly-6A/E (Sca-1) Monoclonal Antibody For Fc - P05533

Item Number: CM0008951

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse Ly-6A/E (Sca-1), a GPI-anchored T-cell activating protein. Suitable for flow cytometry applications. Supports stem cell and immunology research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse Ly-6A/E (Sca-1) mAb
Remarks/Alias	TAP; Sca1; Sca-1; Ly-6A.2; Ly-6A/E; Ly-6E.1
Species	Mouse
GeneID	110454
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-111 of mouse Ly-6A/E (Sca-1) (NP_034868.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P05533
Protein Weight	114kDa
Shipping	Ice bag