

Pe Rabbit Anti-Mouse Ly-6G Monoclonal Antibody For Fc - P35461

Item Number: CM0008787

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse Ly-6G (Gr-1), suitable for flow cytometric detection of mouse granulocytes. GPI-anchored bone marrow antigen, ideal for immunophenotyping and diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse Ly-6G mAb
Remarks/Alias	Gr1; Gr-1; Ly-6G
Species	Mouse
GeneID (Human)	80740
GeneID	546644
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-134 of mouse Ly-6G (NP_001297367.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P35461
Protein Weight	28kDa
Shipping	Ice bag