

Pe Rabbit Anti-Mouse Pd-L2/Cd273/B7-Dc Monoclonal Antibody For Flow Cytometry - Q9Wu15

Item Number: CM0009764



Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse PD-L2/CD273/B7-DC, validated for flow cytometry. Functions as a PD-1 ligand regulating T-cell tolerance. Suitable for immunology research.

[Learn More](#)

Parameter	Specification
Product Name	PE Rabbit anti-Mouse PD-L2/CD273/B7-DC mAb
Remarks/Alias	Btdc; B7-DC; PD-L2; F730015O22Rik
Species	Mouse
Gene ID (Human)	80380
Gene ID	58205
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-123 of mouse PD-L2/CD273/B7-DC (NP_067371.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9WUL5
Protein Weight	28kDa
Shipping	Ice bag