

Pe Rabbit Anti-Mouse Rank/Cd265 Monoclonal Antibody For Flow Cytometry - O35305

Item Number: CM0004566



Introduction

PE Rabbit anti-Mouse RANK/CD265 monoclonal antibody for flow cytometry. Recognizes mouse RANK (ODFR), a key receptor in RANKL-mediated osteoclastogenesis and immune regulation. Suitable for bone and immunology studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse RANK/CD265 mAb
Remarks/Alias	OFE; ODFR; Rank; Ly109; mRANK; TRANCE-R
Species	Mouse
GeneID (human)	8792
GeneID	21934
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-213 of mouse RANK/CD265 (NP_033425.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	O35305
Protein Weight	67kDa
Shipping	Ice bag