


Pe Rabbit Anti-Mouse S100A9 Monoclonal Antibody For Intracellular Flow Cytometry - P31725

Item Number: CM0005202



Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse S100A9 (MRP14), validated for intracellular flow cytometry. S100A9 is a calcium-binding protein involved in inflammation and innate immunity. Cross-reacts with mouse. Ideal for immunological research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse S100A9 mAb
Remarks/Alias	p14; Cagb; GAGB; L1Ag; BEE22; MRP14; 60B8Ag
Species	Mouse
GeneID (Human)	6280
GeneID (Mouse)	20202
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of mouse S100a9 (NP_033140.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P31725
Protein Weight	13kDa
Shipping	Ice bag