


Pe Rabbit Anti-Mouse T-Bet/Tbx21 Monoclonal Antibody For Fc (Intra) - Q9Jkd8

Item Number: CM0005838



Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse T-bet/Tbx21, suitable for intracellular flow cytometry (FC). Ideal for studying Th1 lineage commitment, transcriptional regulation, and antiviral responses in mouse models. - Q9JKD8

[Learn More](#)

Product Name	PE Rabbit anti-Mouse T-bet/Tbx21 mAb
Remarks/Alias	TBT1; Tbet; Tblym
Species	Mouse
GeneID (Human)	30009
GeneID	57765
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 321-530 of mouse T-bet/Tbx21 (NP_062380.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	Q9JKD8
Protein Weight	58 kDa
Shipping	Ice bag