

Pe Rabbit Anti-Mouse Tnf-A Monoclonal Antibody For Flow Cytometry (Intracellular) - P06804

Item Number: CM0004824



Introduction

PE-conjugated rabbit monoclonal antibody against mouse TNF- α , validated for intracellular flow cytometry. Detects the pro-inflammatory cytokine TNF- α in macrophages and other immune cells.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse TNF- α mAb
Remarks/Alias	DIF; Tnfa; TNF-a; TNFSF2; Tnlg1f; Tnfsf1a; TNFalpha; TNF-alpha
Species	Mouse
Gene ID (Human)	7124
Gene ID (Mouse)	21926
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 80-235 of mouse TNF- α (NP_038721.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P06804
Protein Weight	26kDa
Shipping	Ice bag