

Pe Rabbit Anti-Mouse/Rat Cd27 Monoclonal Antibody For Fc - P41272

Item Number: CM0012438

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse/rat CD27 (TNFRSF7). Validated for flow cytometry, reacts with mouse and rat. Suitable for T-cell, B-cell, and NK-cell immunophenotyping studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse/Rat CD27 mAb
Remarks/Alias	S152; Tp55; Tnfrsf7
Species	Mouse
GenelD (Human)	939
GenelD	21940
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse, Rat
SWISS	P41272
Protein Weight	28kDa
Shipping	Ice bag