

Pe Rabbit Anti-Mouse/Rat Cd31/Pecam1 Mab - Q08481

Item Number: CM0010322

Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse and rat CD31/PECAM1, validated for flow cytometry. Ideal for studying leukocyte transendothelial migration, endothelial cell adhesion, and atherosclerosis research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse/Rat CD31/PECAM1 mAb
Remarks/Alias	Cd31; Pecam; PECAM-1
Species	Mouse
GeneID (human)	5175
GeneID	18613
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-220 of mouse CD31 (NP_001292086.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse, Rat
SWISS	Q08481
Protein Weight	69kDa/77kDa/80kDa/81kDa
Shipping	Ice bag