

## Pe Rabbit Anti-Mouse/Rat Cd54 Mab - P13597

Item Number: CM0012334

### Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse/rat CD54 (ICAM-1, Ly-47) for flow cytometry. Epitope within amino acids 28-485. Cross-reacts with mouse and rat. Suitable for leukocyte adhesion and immune cell migration research.

[Learn More](#)

| Parameter        | Value                                                                                                                               |
|------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PE Rabbit anti-Mouse/Rat CD54 mAb                                                                                                   |
| Remarks/Alias    | CD54; Ly-47; Icam-1; MALA-2                                                                                                         |
| Species          | Mouse                                                                                                                               |
| Gene ID (Human)  | 3383                                                                                                                                |
| Gene ID          | 15894                                                                                                                               |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 28-485 of mouse CD54(NP_034623.1) |
| Source           | Rabbit                                                                                                                              |
| Category         | Monoclonal Antibodies                                                                                                               |
| Application      | FC                                                                                                                                  |
| Cross Reactivity | Mouse, Rat                                                                                                                          |
| SWISS            | P13597                                                                                                                              |
| Protein Weight   | 59kDa                                                                                                                               |
| Shipping         | Ice bag                                                                                                                             |