

Pe Rabbit Anti-Pig Cd45 Monoclonal Antibody (Fc)

Item Number: CM0011087

Introduction

PE-conjugated rabbit monoclonal antibody targeting pig CD45. With an immunogen sequence covering amino acids 27-547, this antibody is validated for flow cytometry (FC) and is ideal for pig immune cell analysis and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Pig CD45 mAb
Remarks/Alias	CD45
Species	Pig
Gene ID (Human)	5788
Gene ID (Pig)	100522631
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-547 of pig CD45 (XM_003130596.6)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Pig
Shipping	Ice bag