

Pe Rabbit Anti-Rat Cd25 Monoclonal Antibody For Flow Cytometry - P26897

Item Number: CM0010883

Introduction

PE-conjugated rabbit monoclonal antibody targeting rat CD25 (IL-2R α , P26897) for flow cytometry. Suitable for Treg and immune tolerance research, derived from recombinant extracellular domain immunogen.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Rat CD25 mAb
Alias	IL2RAC
Species	Rat
Human Gene ID	3559
Gene ID	25704
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 22-235 of rat CD25 (NP_037295.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Rat
SWISS Accession	P26897
Protein Weight	31 kDa
Shipping	Ice bag