

Pe Rabbit Anti-Rat Cd8A Monoclonal Antibody For Fc - P07725

Item Number: CM0012535

Introduction

PE Rabbit anti-Rat CD8a monoclonal antibody for flow cytometry (FC). Specific to rat CD8a (32 kDa), a T-cell coreceptor essential for MHC class I recognition and immune response. Ideal for immunology research and diagnostics development.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Rat CD8a mAb
Remarks/Alias	T8; Lyt2; Ly-2; CD8α
Species	Rat
GeneID (Human)	925
GeneID	24930
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Rat
SWISS	P07725
Protein Weight	32 kDa
Shipping	Ice bag