

Pe Rabbit Anti-Rat Ifn-Gamma Monoclonal Antibody For Intracellular Flow Cytometry - P01581

Item Number: CM0010879

Introduction

PE Rabbit anti-Rat IFN-gamma monoclonal antibody validated for intracellular flow cytometry. Recognizes the mature 18 kDa interferon gamma protein. Immunogen: recombinant fragment (aa 20-156). No cross-reactivity with other species. Suitable for rat immune cell analysis.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Rat IFN-gamma mAb
Remarks/Alias	If2f, IFNG2, IFN gamma
Species	Rat
Gene ID (Human)	3458
Gene ID (Rat)	25712
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-156 of rat IFN gamma (NP_620235.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Rat
SWISS	P01581
Protein Weight	18 kDa
Shipping	Ice bag