

Pe Rabbit Igg Isotype Control - Flow Cytometry

Item Number: CM0010576

Introduction

PE-conjugated rabbit IgG isotype control monoclonal antibody for flow cytometry (FC) applications. Designed to assess non-specific binding in rabbit-derived antibody experiments. Ships on ice.

[Learn More](#)

Parameter	Specification
Product Name	PE Rabbit IgG isotype control
Immunogen	Synthetic peptide derived from rabbit IgG isotype control
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Shipping	Ice bag