

Pe Rat Anti-Mouse Cd3 Monoclonal Antibody For Flow Cytometry - P04235 / P22646

Item Number: CM0008507

Introduction

PE-conjugated rat anti-mouse CD3 monoclonal antibody for flow cytometry. Recognizes mouse CD3 complex (delta/epsilon chains). Ideal for T-cell immunophenotyping in murine immune studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rat anti-Mouse CD3 mAb
Remarks/Alias	T3d; CD3; T3e; CD3epsilon
Species	Mouse
Gene ID (Human)	916/915
Gene ID	12500/12501
Immunogen	A synthesized peptide derived from mouse CD3
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P04235/P22646
Protein Weight	19 kDa/21 kDa
Shipping	Ice bag