

Pe Rat Anti-Mouse Cd45R/B220 Monoclonal Antibody For Flow Cytometry - P06800 / P08575

Item Number: CM0006027

Introduction

PE-conjugated rat monoclonal antibody targeting mouse CD45R/B220 (leukocyte common antigen). Validated for flow cytometry, cross-reacts with mouse samples. Ideal for immune cell phenotyping in research and diagnostic development.

[Learn More](#)

Parameter	Description
Product Name	PE Rat anti-Mouse CD45R/B220 mAb
Remarks/Alias	Ioc; B220; Cd45; L-CA; Ly-5; T200; CD45R; Lyt-4
Species	Mouse
GeneID (Human)	5788
GeneID	19264/5788
Immunogen	A synthesized peptide derived from mouse/human CD45R/B220
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06800/P08575
Protein Weight	130kDa/140kDa/144kDa
Shipping	Ice bag