

Pe Rat Anti-Mouse Ly-6G Mab For Fc - P35461

Item Number: CM0008426

Introduction

PE-conjugated rat anti-mouse Ly-6G monoclonal antibody for flow cytometry (FC). Targets Lymphocyte antigen 6G, cross-reacts with mouse, molecular weight 14 kDa. Ideal for bone marrow immunophenotyping research. UniProt P35461.

[Learn More](#)

Parameter	Value
Product Name	PE Rat anti-Mouse Ly-6G mAb
Remarks/Alias	Gr1; Gr-1; Ly-6G Ly-6G.1
Species	Mouse
GeneID (human)	80740
GeneID	546644
Immunogen	A synthesized peptide derived from mouse Ly-6G
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P35461
Protein Weight	14kDa
Shipping	Ice bag