

Pe Rat Anti-Mouse Ly-6G/Ly-6C (Gr-1) Monoclonal Antibody For Fc - P35461 / P0Cw02

Item Number: CM0011108

Introduction

PE-conjugated rat anti-mouse Ly-6G/Ly-6C (Gr-1) monoclonal antibody for flow cytometry of mouse granulocytes. Targets bone marrow-expressed GPI-anchored membrane proteins. Validated in mouse samples.

[Learn More](#)

Parameter	Value
Product Name	PE Rat anti-Mouse Ly-6G/Ly-6C (Gr-1) mAb
Remarks/Alias	Gr1; Gr-1; Ly-6G
Species	Mouse
GenelD (Human)	8600
GenelD	546644/17067
Immunogen	Recombinant protein of mouse Ly-6G/Ly-6C (Gr-1)
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P35461/P0CW02
Protein Weight	14kDa
Shipping	Ice bag