

Pe Streptavidin Monoclonal Antibody (Rabbit) - Sa V1, Sa V2

Item Number: CM0003974

Introduction

PE-conjugated Streptavidin monoclonal antibody (Rabbit) for flow cytometry applications. Species-independent cross-reactivity supports a wide range of diagnostic assay development, available in bulk from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	PE Streptavidin
Remarks/Alias	SA V1, SA V2
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Species independent
Shipping	Ice bag