

Pea15 Rabbit Polyclonal Antibody For Wb And Elisa - Q15121

Item Number: CM0030532



Introduction

PEA15 Rabbit pAb is a polyclonal antibody raised against recombinant human PEA15 (aa 1-130). Suitable for WB and ELISA, reacts with human and mouse. Use for apoptosis, glucose transport and MAP kinase signaling research.

[Learn More](#)

Parameter	Value
Product Name	PEA15 Rabbit pAb
Remarks/Alias	PED; MAT1; HMAT1; MAT1H; PEA-15; HUMMAT1H; PED-PEA15; PED/PEA15; PEA15
Species	Human
GeneID (Human)	8682
GeneID	8682
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human PEA15 (NP_003759.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15121
Protein Weight	15kDa
Shipping	Ice bag