

Pea3 / Etv4 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P43268

Item Number: CM0027452



Introduction

Pea3 / ETV4 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat ETV4, a transcription factor involved in neurogenesis and keratinocyte differentiation. Ideal for cancer and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	Pea3 / ETV4 Rabbit pAb
Remarks/Alias	E1AF; PEA3; E1A-F; PEAS3; Pea3 / ETV4
Species	Human
GeneID (Human)	2118
GeneID	2118
Immunogen	Recombinant Protein corresponding to a sequence within amino acids 1-284 of human ETV4 (NP_001977.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43268
Protein Weight	54kDa
Shipping	Ice bag