

PE/Cyanine5 Mouse Anti-Human Cd3 Monoclonal Antibody For Flow Cytometry - P07766

Item Number: CM0004444



Introduction

PE/Cyanine5-conjugated anti-human CD3 epsilon monoclonal antibody validated for flow cytometry. Suitable for T-cell identification and immune response studies in human samples.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine5 Mouse anti-Human CD3 mAb
Remarks/Alias	T3E, TCRE, IMD18, CD3epsilon
Species	Human
GeneID (human)	916
Immunogen	Recombinant protein of human CD3
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag