


PE/Cyanine5 Rabbit Anti-Mouse Cd11B Monoclonal Antibody For Flow Cytometry - P05555

Item Number: CM0005252



Introduction

High-performance PE/Cyanine5-conjugated rabbit anti-mouse CD11b monoclonal antibody validated for flow cytometry. Targets mouse integrin alpha-M (Itgam), a receptor for complement iC3b, fibrinogen, and ICAM1, crucial for immune cell adhesion and phagocytosis studies.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine5 Rabbit anti-Mouse CD11b mAb
Remarks/Alias	CR3; CR3A; MAC1; Cd11b; Ly-40; Mac-1; Mac-1a; CD11b/CD18; F730045J24Rik
Species	Mouse
GeneID (Human)	3684
GeneID (Mouse)	16409
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 143-337 of mouse CD11b (NP_032427.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P05555
Protein Weight	127kDa
Shipping	Ice bag