

Pe/Cyanine5 Rabbit Anti-Mouse Cd5 Monoclonal Antibody For Flow Cytometry - P13379

Item Number: CM0007872

Introduction

PE/Cyanine5-conjugated rabbit monoclonal antibody targeting mouse CD5 (Ly-1), validated for flow cytometry. Recognizes a 54 kDa transmembrane glycoprotein on T-cells and B1a cells; specific for mouse specimens. Ideal for immunophenotyping and immunological research.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine5 Rabbit anti-Mouse CD5 mAb
Remarks/Alias	Ly-1; Ly-A; Ly-12; Lyt-1
Species	Mouse
GeneID (Human)	921
GeneID (Mouse)	12507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-371 of mouse CD5 (NP_031676.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P13379
Protein Weight	54kDa
Shipping	Ice bag