

## Pe/Cyanine5 Rabbit Anti-Mouse Cd90.2 Mab For Flow Cytometry - P01831

Item Number: CM0010248

### Introduction

PE/Cyanine5-conjugated rabbit monoclonal antibody specific for mouse CD90.2 (Thy-1), suitable for flow cytometry. Detects a GPI-anchored membrane glycoprotein implicated in cell-cell interactions during synaptogenesis. Ideal for immunological studies.

[Learn More](#)

| Parameter        | Value                                                                |
|------------------|----------------------------------------------------------------------|
| Product Name     | PE/Cyanine5 Rabbit anti-Mouse CD90.2 mAb                             |
| Remarks/Alias    | T25; CD90; Thy-1; Thy1.1; Thy1.2; Thy-1.2                            |
| Species          | Mouse                                                                |
| GeneID (human)   | 7070                                                                 |
| GeneID           | 21838                                                                |
| Immunogen        | Synthetic peptide   A synthesized peptide derived from mouse CD90.2. |
| Source           | Rabbit                                                               |
| Category         | Monoclonal Antibodies                                                |
| Application      | FC                                                                   |
| Cross Reactivity | Mouse                                                                |
| SWISS            | P01831                                                               |
| Protein Weight   | 18kDa                                                                |
| Shipping         | Ice bag                                                              |