

PE/Cyanine7 Anti-Human Cd27 Monoclonal Antibody For Flow Cytometry - P26842

Item Number: CM0008504



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody against human CD27 (TNFRSF7), a costimulatory immune checkpoint receptor. Validated for flow cytometry (FC) applications. Ideal for T-cell and immune profiling research.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human CD27 mAb
Remarks/Alias	T14; S152; Tp55; TNFRSF7; S152. LPFS2
Species	Human
GeneID (human)	939
GeneID	939
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-192 of human CD27 (NP_001233.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P26842
Protein Weight	29kDa
Shipping	Ice bag