

PE/Cyanine7 Anti-Mouse Cd11C Monoclonal Antibody For Fc - Q9Qxh4

Item Number: CM0007893

Introduction

PE/Cyanine7-conjugated monoclonal antibody targeting mouse CD11c (Integrin alpha-X, Q9QXH4), validated for flow cytometry. Recognizes fibrinogen receptor involved in inflammatory cell adhesion and monocyte chemotaxis.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 anti-Mouse CD11c mAb
Remarks/Alias	Cr4; N418; Cd11c
Species	Mouse
GeneID (Human)	3687
GeneID	16411
Immunogen	Recombinant protein of mouse CD11c
Source	Armenian Hamster
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9QXH4
Protein Weight	129kDa
Shipping	Ice bag