

PE/Cyanine7 Anti-Mouse Cd161C/Nk1.1 Monoclonal Antibody For Flow Cytometry - P27814

Item Number: CM0008238

Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody specifically recognizing mouse CD161c/NK1.1 (NKR-P1C). Validated for flow cytometry (FC) to study NK cell function and immunophenotyping. Suitable for activation and cytotoxicity assays.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Mouse CD161c/NK1.1 mAb
Remarks/Alias	Nk1, Ly59, Nk-1, CD161, Ly-59, Ly55c, NK1.1, NKR-P1, Nk1.2, Klrb1b, NK-RP1, Nk-1.2, Nkrp1b, Nkrp1c, ly-55c, NKR-P1.9
Species	Mouse
GeneID (Human)	3820
GeneID	17059
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 67-223 of mouse CD161/NK1.1 (NP_001153376.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P27814
Protein Weight	25kDa
Shipping	Ice bag