

Pe/Cyanine7 Anti-Mouse Cd24 Monoclonal Antibody For Flow Cytometry - P24807

Item Number: CM0010755

Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting mouse CD24, optimized for flow cytometry. Use in research on cell differentiation, thymocyte development, B-cell activation, and autoimmunity.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Mouse CD24 mAb
Remarks/Alias	HSA; Cd24; Ly-52; nectadrin
Species	Mouse
GenelD (Human)	100133941
GenelD	12484
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-53 of mouse CD24 (NP_033976.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P24807
Protein Weight	8kDa
Shipping	Ice bag