

PE/Cyanine7 Anti-Mouse Cd3E Monoclonal Antibody For Flow Cytometry - P22646

Item Number: CM0009743



Introduction

PE/Cyanine7-conjugated Armenian Hamster anti-Mouse CD3e monoclonal antibody for flow cytometry detection of mouse CD3e, a T-cell surface glycoprotein involved in adaptive immunity - P22646.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Armenian Hamster anti-Mouse CD3e mAb
Remarks/Alias	CD3; T3e; CD3epsilon
Species	Mouse
GeneID (Human)	916
GeneID	12501
Immunogen	Recombinant protein
Source	Armenian Hamster
Category	Monoclonal Antibodies
Application	Flow Cytometry
Cross Reactivity	Mouse
SWISS	P22646
Protein Weight	21kDa
Shipping	Ice bag