

PE/Cyanine7 Anti-Mouse Ly-6A/E (Sca-1) Monoclonal Antibody For Flow Cytometry - P05533

Item Number: CM0004378

Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting mouse Ly-6A/E (Sca-1) for flow cytometry. Validated for mouse reactivity; ideal for stem cell and T-cell activation studies. One-stop sourcing from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Mouse Ly-6A/E (Sca-1) mAb
Remarks/Alias	TAP; Sca1; Sca-1; Ly-6A.2; Ly-6A/E; Ly-6E.1
Species	Mouse
Gene ID	110454
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-111 of mouse Ly-6A/E (Sca-1) (NP_034868.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P05533
Protein Weight	115kDa
Shipping	Ice bag