

Pe/Cyanine7 Anti-Mouse Tcr B Chain Mab For Fc - A0A0B4J1G9

Item Number: CM0011133

Introduction

PE/Cyanine7-conjugated monoclonal antibody against mouse TCR β chain (Armenian Hamster, FC-validated). Reliable for flow cytometry detection of TCR β in mouse samples. - A0A0B4J1G9

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 anti-Mouse TCR β chain mAb
Remarks/Alias	Tib; TCRbeta
Species	Mouse
GeneID (human)	28614
GeneID	21577
Immunogen	Synthetic peptide A synthesized peptide derived from mouse TCR β chain.
Source	Armenian Hamster
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	A0A0B4J1G9
Protein Weight	13kDa
Shipping	Ice bag