

PE/Cyanine7 Mouse Anti-Human Cd11B Monoclonal Antibody For Flow Cytometry - P11215

Item Number: CM0012509



Introduction

PE/Cyanine7-conjugated mouse monoclonal antibody targeting human CD11b (Integrin alpha-M). Validated for flow cytometry, this antibody binds specifically to a 127 kDa myeloid receptor expressed on monocytes and granulocytes. Ideal for inflammation and immunology research applications.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Mouse anti-Human CD11b mAb
Remarks/Alias	CR3A; MO1A; CD11B; MAC-1; MAC1A; SLEB6
Species	Human
Gene ID	3684
Immunogen	A synthesized peptide derived from human CD11b
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P11215
Protein Weight	127kDa
Shipping	Ice bag