

PE/Cyanine7 Mouse Anti-Human CD11C Monoclonal Antibody For Fc - P20702

Item Number: CM0009887



Introduction

PE/Cyanine7-conjugated mouse monoclonal antibody specific for human CD11c (Integrin alpha-X, 128 kDa). Suitable for flow cytometry, recognizes a fibrinogen receptor on monocytes/granulocytes involved in inflammation and erythrophagocytosis. - P20702

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Mouse anti-Human CD11c mAb
Remarks/Alias	CD11C; SLEB6
Species	Human
GeneID	3687 (Human)
Immunogen	Recombinant protein of human CD11c
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P20702
Protein Weight	128kDa
Shipping	Ice bag