

Pe/Cyanine7 Mouse Anti-Human Cd45Ro Mab - P08575

Item Number: CM0008459



Introduction

High-quality PE/Cyanine7 Mouse anti-Human CD45RO monoclonal antibody for flow cytometry (FC). Recognizes human CD45RO isoform of PTPRC, a key T-cell activation marker. Ideal for immunological research and clinical diagnostics.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Mouse anti-Human CD45RO mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
GeneID (Human)	5788
GeneID	5788
Immunogen	Synthetic peptide derived from Human CD45RO
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08575
Protein Weight	136kDa
Shipping	Ice bag