


PE/Cyanine7 Rabbit Anti-Human Axl Monoclonal Antibody For Flow Cytometry - P30530

Item Number: CM0007924



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody against human AXL (UFO). Suitable for flow cytometry (FC) on human samples. AXL is a receptor tyrosine kinase involved in cell survival, proliferation, and virus entry. Protein weight: 98 kDa.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human AXL mAb
Remarks/Alias	ARK; UFO; JTK11; Tyro7
Species	Human
GeneID (Human)	558
GeneID	558
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-449 of human AXL (NP_068713.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P30530
Protein Weight	98kDa
Shipping	Ice bag