

PE/Cyanine7 Rabbit Anti-Human Cd13/Anpep Mab For Fc - P15144

Item Number: CM0011856



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting human CD13/ANPEP, validated for flow cytometry. Ideal for immunophenotyping, tumor vasculature analysis, and viral receptor studies. - P15144

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human CD13/ANPEP mAb
Remarks/Alias	APN; AP-M; AP-N; CD13; LAP1; P150; PEPN; hAPN; GP150
Species	Human
Gene ID (Human)	290
Gene ID	290
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 69-967 of human CD13/ANPEP (NP_001141.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P15144
Protein Weight	110kDa
Shipping	Ice bag