


Pe/Cyanine7 Rabbit Anti-Human Cd2 Monoclonal Antibody For Fc - P06729

Item Number: CM0006000



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody against human CD2 (T11/LFA-2) for flow cytometry, validated for human samples. Ideal for T-cell and NK cell adhesion and activation studies.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human CD2 mAb
Remarks/Alias	T11; SRBC; LFA-2
Species	Human
GeneID (human)	914
GeneID	914
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-209 of human CD2 (NP_001758.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06729
Protein Weight	39kDa
Shipping	Ice bag