


PE/Cyanine7 Rabbit Anti-Human Cd314 Monoclonal Antibody For Fc - P26718

Item Number: CM0005389



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting human CD314 (NKG2D) for flow cytometry. Recognizes a 25kDa type II membrane protein involved in NK cell-mediated immunosurveillance. Suitable for human samples.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human CD314 mAb
Remarks/Alias	KLR; CD314; NKG2D; NKG2-D; D12S2489E
Species	Human
GeneID (Human)	22914
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 73-216 of human CD314 (NP_031386.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P26718
Protein Weight	25kDa
Shipping	Ice bag