


PE/Cyanine7 Rabbit Anti-Human Cd33 Monoclonal Antibody For Flow Cytometry - P20138

Item Number: CM0004511



Introduction

PE/Cyanine7-conjugated rabbit anti-human CD33 monoclonal antibody for flow cytometry. Recognizes CD33 (Siglec-3, p67) on monocytic/myeloid cells. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human CD33 mAb
Remarks/Alias	p67; SIGLEC3; SIGLEC-3
Species	Human
GeneID (Human)	945
GeneID	945
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-259 of human CD33(NP_001763.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P20138
Protein Weight	40kDa
Shipping	Ice bag