


Pe/Cyanine7 Rabbit Anti-Human Cd38 Mab For Fc - P28907

Item Number: CM0008346



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody against human CD38 for flow cytometry. CD38 is a multifunctional ectoenzyme and receptor involved in calcium signaling and lymphocyte migration. Species reactivity: Human. SWISS: P28907.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human CD38 mAb
Remarks/Alias	ADPRC1; cADPR1; ADPRC 1
Species	Human
GeneID (human)	952
GeneID	952
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 43-300 of human CD38 (P28907).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P28907
Protein Weight	13kDa/34kDa
Shipping	Ice bag