

## Pe/Cyanine7 Rabbit Anti-Human Cd42B Monoclonal Antibody For Fc - P07359

Item Number: CM0008377



### Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting human CD42b (platelet glycoprotein Ib alpha chain) for flow cytometry. Recognizes native GP1BA protein on human platelets, ideal for platelet adhesion and hemostasis research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | PE/Cyanine7 Rabbit anti-Human CD42b mAb                                                                           |
| Remarks/Alias    | BSS, GP1B, VWDP, CD42B, GPIbA, BDPLT1, BDPLT3, DBPLT3, GPIbalpha, CD42b-alpha                                     |
| Species          | Human                                                                                                             |
| Gene ID          | 2811                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 15-531 of human CD42b (NP_000164.5) |
| Source           | Rabbit                                                                                                            |
| Category         | Monoclonal Antibodies                                                                                             |
| Application      | FC                                                                                                                |
| Cross Reactivity | Human                                                                                                             |
| SWISS            | P07359                                                                                                            |
| Protein Weight   | 72 kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                           |