

Pe/Cyanine7 Rabbit Anti-Human Cd45Rb Monoclonal Antibody For Flow Cytometry - P08575-9

Item Number: CM0012415

Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting human CD45RB (isoform 7 of PTPRC). Recommended for flow cytometry. Detects leukocyte common antigen, suitable for immunophenotyping of human samples.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human CD45RB mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
Gene ID (Human)	5788
Gene ID	5788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-463 of human CD45RB (XP_006711537.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08575-9
Protein Weight	135kDa
Shipping	Ice bag