


PE/Cyanine7 Rabbit Anti-Human Cd5 Monoclonal Antibody For Flow Cytometry - P06127

Item Number: CM0005499



Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting human CD5 (T1/Leu-1), a pan-T-cell and B1a-cell surface receptor. Validated for flow cytometry, ideal for immunophenotyping and T-cell signaling studies.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human CD5 mAb
Remarks/Alias	T1; LEU1
Species	Human
GeneID	921
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-371 of human CD5 (NP_055022.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06127
Protein Weight	55kDa
Shipping	Ice bag