

PE/Cyanine7 Rabbit Anti-Human Cd8B Monoclonal Antibody For Fc - P10966

Item Number: CM0010212

Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody targeting human CD8b (P10966) for flow cytometry. Suitable for studying immune response and T-cell coreceptor functions. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human CD8b mAb
Remarks/Alias	LY3; P37; LEU2; LYT3; Ly-3; CD8B1; CD8beta
Species	Human
GeneID (Human)	926
GeneID	926
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-170 of human CD8b (NP_004922.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P10966
Protein Weight	24kDa
Shipping	Ice bag