

PE/Cyanine7 Rabbit Anti-Human Cd90/Thy1 Monoclonal Antibody For Flow Cytometry - P04216

Item Number: CM0004421

Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody against human CD90/Thy1, designed for flow cytometry applications. Ideal for immunophenotyping and research in stem cell biology, neuroscience, and cell signaling.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human CD90/Thy1 mAb
Remarks/Alias	CD90; CDw90
Species	Human
GeneID (Human)	7070
GeneID	7070
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 20-130 of human CD90/Thy1 (NP_006279.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P04216
Protein Weight	18kDa
Shipping	Ice bag