

PE/Cyanine7 Rabbit Anti-Human Cr2/Cd21 Monoclonal Antibody For Flow Cytometry - P20023

Item Number: CM0005392



Introduction

PE/Cyanine7-conjugated rabbit anti-human CR2/CD21 monoclonal antibody validated for flow cytometry. Recognizes the B-cell coreceptor CR2 (CD21, C3d receptor) involved in complement pathways and EBV infection. Suitable for human immunology research.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human CR2/CD21 mAb
Remarks/Alias	CR; C3DR; CD21; CVID7; SLEB9
Species	Human
Gene ID (Human)	1380
Gene ID	1380
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-971 of human CR2/CD21 (NP_001868.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P20023
Protein Weight	110kDa/113kDa/118kDa/119kDa
Shipping	Ice bag